IKA DZM control.m
IKA DZM control.o

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OPERATING INSTRUCTIONS GB 10

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Garantie

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards. In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

Correct use

The IKA DZM-M is in connection with the sensors DZM-S.o and DZM-S.m a measuring system for contactless speed measurement.

Associated standards and regulations

Associated EU guidelines: EMC guidelines 89/336/EEC

Construction in accordance with the following safety standards:

EN 61 010-1 / VDE 411-1	EN 50 081	EN 50 082	EN 55 011
EN 60 055	CAN/CSA C22.2 (1010-1)	UL 3101-1	

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Contents

	Page
Guarantee	3
Correct use	3
Associated standards and regulations	3
Unpacking	10
Safety instructions	10
Setting up the instrument	10
General information	10
Assembly of components	11
Sensor assembly	11
Possible uses	12
Interfaces	13
Configuration of the serial RS232 C Interface	13
Instruction Syntax and Format	14
NAMUR - Instructions	14
PC 5.1 Adapter IKA-Control	14
PC 2.1 PC-Cable	14
Maintenance and cleaning	15
Accessories	15
Technical data	15
List of spare parts	22
Spare parts diagram	23

Unpacking

Unpack the machine carefully and check it thoroughly for damage. It is important that any damage resulting from transit is determined when unpacking. An immediate damage report should be filed if necessary (post office, rail authority or shipping company).

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Safety instructions

When setting up: First connect the plug of the power supply unit to the jack of the DZM-M, then insert the mains plug of the power supply unit into the wall outlet.

When dismantling: First pull the mains plug of the power supply unit from the wall outlet, then disconnect the plug from the jack of the DZM-M if necessary.

Setting up in the instrument

Please pay attention to the ambient conditions listed in Technical data (temperature and humidity).

General information

The IKA DZM-control displays rotational speeds and stroke frequencies in digital form. For speeds > 1999 rpm, the units and tens columns are shifted to the right and displayed on the smaller 7-segment LCD fields.

Because of the measurement principle employed, there may be some delay before the measurement value is correctly displayed (e.g. after a sudden speed increase from 0 to 100 rpm).

To facilitate checking of the test set-up, the green LED is used to provide a pulse display.

The measured values can be recorded via the built-in RS232 interface, the analog output or the frequency output.

The DZM-control can be used in combination with a PC and IKASOFT dc process monitoring and documentation to ISO 9000.

Various speed sensors can be connected to the DZM-M monitor:

DZM-S.o: Optical sensor (photoelectric barrier principle)

DZM-S.m: Magnetic sensor (Hall principle)

Assembly of components

The monitor for reading the indicated speeds can be used as either a bench-top or a stand-mounted unit. First connect the power supply unit to the monitor by plugging the jack plug into the jack socket provided on the monitor. Then connect a sensor to the monitor (further details below).

Connect any necessary adaptor cable, RS232 interface cable to the PC, analog output cable etc. Only now should the power supply unit be plugged into the socket. The measuring system is switched on with the rocker switch on the left side of the monitor housing.

Assembly of sensor

If using the DZM-S.o sensor, ensure that no stray light or water spray can enter the sensor tube.

Sensor DZM-S.o

Mounting on IKA agitators or IKA dispersion units with sensor socket:

In this case simply remove the protective cap from the sensor socket and insert the measuring sensor into the socket.

Mounting on an agitator without sensor socket:

In this case optical markings (e.g. the matt black adhesive tape included) must be affixed to the chuck of the agitator.

The DZM-control can recognize and process either 1 or 10 markings around the circumference.

If the specified dimensions are maintained, the DZM-control automatically detects whether 1 or 10 markings are affixed around the circumference and accordingly displays the correct rotational speed with a resolution of 1 rpm.

With 10 markings, the light / dark scanning ratio must be $> 1/3$ and $< 2/3$, ideally 1/1.

For speeds < 2000 rpm it is advisable to affix 10 markings as evenly spaced as possible around the circumference.

Example 1

You have a drill chuck with a diameter of 34mm and wish to affix 10 markings.

Choosing the width of the markings:

Circumference:	$U = \pi \times d$	$U = 107\text{mm}$
Part width:	$bt = U / 10$	$bt = 10,7$
Marking width:	$bm^{\min} = 1/3 \times bt$	$bm^{\min} = 3,6\text{mm}$
	$bm^{\max} = 2/3 \times bt$	$bm^{\min} = 7,1\text{mm}$

Selected marking width $bm = 5\text{mm}$.

With 1 marking, the light / dark scanning ratio must be $< 1/5$ or $> 4/5$.

If only one marking is to be affixed, the light / dark scanning ratio must be $< 1/5$ for speeds < 150 rpm, otherwise the measuring system will be more sensitive to interference.

Example 2

You have a drill chuck with a diameter of 34mm and wish to affix 1 marking.

Choosing the width of the marking:

Circumference:	$U = \pi \times d$	$U = 107\text{mm}$
Part width:	$bt = U / 1$	$bt = 107$
Marking width:	$bm^{\min} = 4/5 \times bt$	$bm^{\min} = 85,6\text{mm}$
	$bm^{\max} = 1/5 \times bt$	$bm^{\min} = 21,4\text{mm}$

Selected marking width for speeds < 150 rpm $bm = 90\text{mm}$ and for speeds > 150 rpm $bm = 10\text{mm}$.

The DZM-S.o sensor is fixed, with a lab clamp for example, and centered with respect to the drill chuck. The distance between the drill chuck and the edge of the sensor must be 5-10mm.

Because of the automatic detection of the number of markings, an excessively high speed may be displayed for a short time before the marking number is detected.

Sensor DZM-S.m

Mounting on IKA magnetic agitators

The housing of the DZM-S.m sensor is shaped to allow it to be fastened directly onto the bar of a lab stand (dia. 10mm). The "eye" of the sensor (oval bump on housing) should be centered on the center of the agitator plate and should be approximately on the same level with it.

The sensor housing and the connecting cable must no touch the heating plate. A minimum distance of 10mm must be maintained. (Temperature of heating plate 300°C)

Mounting on an IKA shaker (e.g. KS 125 basic)

The permanent magnet must be mounted flat on the shaking table (rear edge) using double-sided adhesive tape. The DZM-S.m sensor must be aligned with the magnet block (included) such that the „eye“ of the sensor (oval bump on housing) is aimed at the center of the magnet. A maximum distance of approx. 10mm between the sensor and the magnet must be maintained.

The DZM-S.m sensor (lying on this right side) can be fastened to the shaker housing using double-sided adhesive tape or self-adhesive velcro stripes (included with the DZM-S.m).

Mounting for other applications

The DZM-control can be used for a wide variety of measurement tasks in industry and the laboratory if the procedures described above are observed.

Possible uses

Equipment	DZM-S.m	DZM-S.o	Cable for T50
Color Squid	X		
RCT basic	X		
RET basic	X		
RH	X		
Mini MR standard	X		
Mini MR 1 basic	X		
REO	X		
EUROSTAR basic		X	
EUROSTAR power basic		X	
RW 16 basic		X	
RW 20.n		X	
RW 25		X	
RW 28		X	
RW 47		X	
MS 1 / 2	X		
MTS 2 / 4	X		
KS 125 basic	X		
KS / HS 250 basic	X		
T 20 standard		X	
T 25 basic		X	
T 50			X
SD 41		X	
HBR 1 basic	X		
TER 2	X		
IKA rotation evaporator (Marking on steam pipe)		X	


In addition the unit can also be used with older IKA devices and those of other manufacturers.

Interfaces

There is a 15-pin interface jack on the right-hand side of the DZM-M monitor. This interface is used to transfer analog and digital signals. With the adaptor PC 5.1, this 15-pin interface can be split into an analog and frequency output (round 7-pin connector) and a 9-pin RS 232 interface. A data recorder, for example, can be connected to the frequency output.

The 9-pin RS 232 interface can be connected to a PC. The rotation speed can then be monitored and documented to ISO 9000 with the aid of the IKA software „IKASOFT dc“.

Connecting for measurement sensor, analog output and serial RS 232 interface.

Configuration of the serial RS 232 C Interface

- The functions of the interface lines between laboratory instrument and automation system are selected from the specified signals of the EIA-standard RS 232 C, corresponding with DIN 66 020 part 1.
- For the electrical properties of the interface lines and for the allotment of the signal status, standard RS 232 C, corresponding with DIN 66 259, part 1 applies.
- Transmission-method: Asynchronous signal transmission in start-stop-operation
- Mode of transmission: Fully Duplex
- Character format: Character imaging acc. to data format DIN 66 022 for start-stop-operation
1 start bit
7 character bits
1 parity bit (even)
1 stop bit
- Transmission speed: 9600 bit / s
- Data flow control: Hardware handshake RTS / CTS
RTS(Pin 7)LOW (positive tension): PC may send
RTS(Pin 7)HIGH (negative tension): PC may not send
CTS(Pin 8)LOW (positive tension): PC recipient
CTS(Pin 8)HIGH (negative tension): PC not recipient
- Access method: Data communication from laboratory instrument to processor is only possible on demand of the processor

Instruction Syntax and Format

Here applies the following:

- The instructions are generally sent from the processor (master) to the laboratory instrument (slave).
- The laboratory instrument exclusively sends on demand of the processor. Even error codes cannot be spontaneously communicated from the laboratory instrument to the processor (automation system).
- The instructions are transmitted in capital letters.
- Instructions and parameters as well as subsequent parameters are separated by at least **one** blank. (Code: hex 0x20)
- Each individual instruction including parameters and data as well as each replay are terminated with CR LF abgeschlossen (Code: hex 0x0D 0x0A) and have a maximum length of 80 characters.
- The decimal separator in a floating point is the point (Code hex 0x2E).

The above statements largely correspond with the recommendations of the NAMUR-Association. (NAMUR-recommendations for the design of electric plug connections for the analog and digital signal transmission to laboratory - MSR individual units. Rev. 1.1).

NAMUR instructions


Overview of the NAMUR-Instructions (which the DZM-M understands)

NAMUR-Instructions	Function
IN_PV_4	Reading the real-value
STATUS_4	Display of status 0: manual operation without fault <0: error code: -83: wrong parity -84: unknown instruction -85: wrong instruction sequence -87: not sufficient storage space
IN_TYPE	Requiring identification of laboratory instrument
IN_NAME	Requiring the designation name

OUT_NAME | output of designation name.
name | (max. 6 digits, default: IKADZM)


PC 5.1 Adapter IKA-Control

An adaptor cable is available as an accessory part. It splits the analog and serial signals. The analog output signals are allotted to a 7-polar coupling according to the NAMUR-recommendation; the serial signals to a 9-polar Sub-D-bushing plug (RS 232 C).

PC 2.1 PC-Cable

For connecting the 9-polar bushing with an PC, a PC-cable is required.

- PC 1.2** IKA adaptor for connecting the PC 1.5 to the PC with a 9-pin SUB-D plug
- PC 1.5** IKA cable for connecting the serial RS232 interface (15-pin HD-SUB-D) on the DZM-M to a PC with 9-pin SUB-D plug
- AK 2.1** IKA cable required to connect data recorders, instruments or similar devices to the analog output (only possible via PC 5.1).

- PC 5.1** IKA adaptor **PC 2.1** IKA Cable
- PC 1.2** IKA adaptor **AK 2.1** IKA Cable
- PC 1.5** IKA Cable
- Meßkabel DZM-K** for T50 to connect to T50 ULTRA-TURRAX 1m für DZM-S.o und DZM-S.m.

Maintenance and cleaning

The instrument is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Accessories

DZM-S.o	Optical sensor
DZM-S.m	Magnetic sensor
R 191	Masuring sensor clamp (for DZM-S.o)

Technical data

Voltage supply:	Power supply unit	12 .. 16V DC
Current consumption:	max. 100mA	at $U_v = 12V$ DC
Perm. ON duration:	100%	
Perm. ambiente temperature:	0 .. 40°C	
	IP 42 nach DIN 40 050	
Perm. relativ humidity:	80%	
Measuring range:	0 - 50 000 rpm (DZM-S.o)	
	0 - 10 000 rpm (DZM-S.m)	
Resolution:	1 rpm	
Display:	LCD	
	Up to 1999 rpm: 3 1/2 digit 11mm high	
	Above 1999 rpm: 3 1/2 digit 11mm high and tens and columns 5mm high (shifted to the right)	
Measurement:	$\leq \pm 1$ rpm $\pm 0,4\%$ (if the instructions and recommendations listed under Assembly of Sensor are observed).	
Sensor connection:	6-pin jack as per DIN 45 322	
Actual speed output:	15-pin HD Sub-D port RS 232	
	Internal resistance of instrument $\geq 10k\Omega$	
	Analog: 1V/ 1000 rpm 0...4000 rpm	
	1V/ 10000 rpm >4000 rpm	
	Frequency: 1 puls per revolution or stroke	
Weight:	0,25 kg (DZM-M)	
Dimensions:	80 x 180 x 75 (W x L x H)	

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List of spare parts

IKA DZM control

Item	Designation
1	Plastite screw 3,5 x 12
2	switch
3	PCB logistics
4	clip
5	curved washer
6	rubber base

8	PCB plug 8/15pin
9	PCB amplifier
10	front panel
12	housing - upper part
13	locking mechanism
14	flat ribbon conductor
15	flat ribbon conductor
16	cover
20	reference plate
2001	housing - upper part, complete

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
IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.
2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works do Brasil Ltda.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

South America

IKA® Works do Brasil Ltda.
Estrada do Guerenguê, 491
Taquara Jacarepagua, RJ
Rio de Janeiro
22713-000 Brasil
TEL. +55 21 2435-9600
FAX +55 21 2435-9601
E-mail: brasil@ika.net

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

IKA® Works Guangzhou
173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn