
Dispergierwerkzeuge/ Dispersion tools/ Outils dispersants

S 25 NK - 19 G

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Safety instructions

- ⚠ - No liquid may get into the drive unit of the disperser.
- ⚠ - The work has to be interrupted immediately if you notice unusual noise and/ or increasing emission of liquid from the upper rinsing drilling. The slide bearing must be controlled, and replaced if necessary.

- While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles. Follow the health and safety regulations accident prevention regulations.
- Note the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.

Product information

- The shaft bearing of the dispersion tool consists of a slide bearing (PTFE) and a ball bearing.
- Content of package: dispersion tool, complete and tool “bar” (see fig. 1)
- The plastic materials used are approved by FDA (Food and Drug Administration Approval).

Application instructions

- Attention! Never run the dispersion tool dry. Without cooling by the medium the bearing will be destroyed.
- The dispersion tool is not suitable for continuous operation.
- Dispersion causes the medium to heat up.
- Likewise the optimal dispersion duration and rotating frequency must be determined by attempts.
- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator (integrated in the shaft tube). Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Designation	Application	Max. permissible number of revolutions [1/min]	Working range [ml]
S 25 NK-19G	Suspension, emulsifying, disagglomeration (batch operation)	24000	25-1500

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Spare parts list

Item	Designation	Spare parts diagram see front page.
1	Shaft tube	
2	Shaft	For orders of spare parts please specify the fabrication number on the shaft tube, type of dispersion tool and the designation of the spare part.
6	Bearing	
7	Rotor	

Cleaning, disinfection and sterilization

Allowed procedures	Methods of sterilisation	Notice
Damp heat	Autoclave to 121 °C at 2 bar positive pressure	Life span of the bearings shortens
Chemical procedures	By germ - killing dissolvents alcohol, phenol,formalin.... Disinfectant remainders are to be removed with germ-free water	Note the chemical resistance of ball bearings ball bearings and PTFE-seal.

For pre-cleaning the dispersion tool, operate it in a dissolvent for dissolving the residuaries of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be disassembled and cleaned immediately after working, so that the adhering residuaries of the substance don't cultivate unwanted bacterial cultures.

For cleaning purposes disassemble the dispersion tool as described below (items see front page):

- Unscrew the rotor (7) from the shaft (2) using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
- Pull the shaft (2) out of the shaft tube in an upward direction. Push the bearing (6) out of the shaft tube (1) by using the bar included with the tool set.
- The bearing (6) shall be replaced in the case of wear.

The assembly of the dispersion tool takes place in reverse order.

Attention! The rotor must not be tightened too strongly, because otherwise the small teeth at the rotor become bent and so the dispersion tool gets destroyed.

Technical Data											
Designation	Working - range ml	Stator ø	Rotor ø	Peripheral speed at 24000 1/min	Min./max. im-mersion depth	Materials in contact with medium	Suitable for solvents	Max. tem-perature range upto	Ultimate fineness dispersions	Ultimate fineness emulsions	
S 25 NK-19 G	25- 1500	19 mm	12,7 mm	15,9 m/s	40/165 mm	PTFE AISI 316L	yes	120 °C	10-50 µm	1-10 µm	

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