
IKA LR 2000 V

| | | |
|-------------------------------|-----------|-----------|
| <i>BETRIEBSANLEITUNG</i> | <i>DE</i> | <i>2</i> |
| <i>OPERATING INSTRUCTIONS</i> | <i>EN</i> | <i>10</i> |
| <i>MODE D'EMPLOI</i> | <i>FR</i> | <i>18</i> |

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| | Page |
|--|------|
| Safety instructions | 10 |
| Intended use | 11 |
| Unpacking | 12 |
| Setting up | 12 |
| Commissioning | 13 |
| Description | 13 |
| Device base, telescopic stand | 13 |
| Stirring motor carrier, stirring shaft bearing | 13 |
| Stirrer drive | 14 |
| Clamping device of vessel | 14 |
| Reactor cover | 14 |
| Use of dispersing tool S 25 KV | 14 |
| Reactor vessel | 15 |
| Materials contacting the product | 15 |
| Cleaning and sterilization | 15 |
| Maintenance | 16 |
| Spare parts | 16 |
| Wear parts | 16 |
| Accessories | 17 |
| Warranty | 17 |
| LR 2000 P | 17 |

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed.
- Work may only be carried out behind a protective screen or under the closed laboratory extractor hood!
- The voltage stated on the nameplate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- There must be an EMERGENCY STOP facility which disconnects all poles to the laboratory reactor from the power supply.
- The EMERGENCY STOP facility must be easy to reach in a danger situation and be outside the area at risk.
- The laboratory reactor must only be operated using the safety switch included in the kit contents.
- The power supply to the agitator and the disperser drive must be through an in-built safety switch. Only ever connect these units through it. That way the agitator and disperser drives can only operate when the vessel tongs are closed. Switch off the agitator and disperser drive before removing the reactor vessel!
- Equipment which have to continue to operate after the reactor lid has been opened (e.g. thermostat) must not be wired through the safety switch! Heating without agitating is not permissible! Warning delay in boiling!
- Working with inflammable or explosive gases, vapours or solvents is not permissible in the laboratory reactor LR2000!
- Do not process any flammable or combustible materials.
- **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Operation with over pressure is not permissible! The reactor vessel may only be operated with a maximum vacuum of 25 mbar.
- Observe correct handling procedures for reaction educts and products.
- Depending on the application and materials, contact with or breathing in poisonous liquids, gases, mists, steam or dust can be dangerous. If relevant, biological or microbiological materials can also be dangerous!
- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact **IKA** if you have any questions.
- Do not touch rotating parts during operation! The anchor agitator and the dispersion tool must not be operated when the reactor vessel is open.
- Starting an open laboratory reactor (even for maintenance purposes) is not permissible!
- Rotating the anchor agitator or dispersal tool with some of the product sticking to it when open can lead to particles or liquid being sprayed.
- The unit may not be operated without a reactor vessel.

The reactor vessels LR 2000.1 and LR2000.2 may only be heated up without pressure through the double jacket with thermostats or comparable heat sources and suitable liquids at constant temperature! Heating with other heat sources is not permissible! Danger of cracking!

- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Protect the appliance and accessories from bumps and impacts.
- When working at normal pressure the reactor system must always be vented so as to prevent an increase in pressure from highly volatile gases or unknown pressure curve during the reaction. Condense volatile gases with a cooler with a conical joint (e.g. return condenser) in the reactor lid.
- Observe the permissible maximum and minimum temperatures of $-10\text{ }^{\circ}\text{C}$ to $+230\text{ }^{\circ}\text{C}$ in the reactor vessel! There is a risk of burning and scalding by the hot reactor vessel or reactor lid! Wear safety gloves when handling heated parts.
- Check that the thermostats for heating are working properly! A defective thermostat could cause the reaction to proceed in an uncontrolled manner.
- Before filling the reactor vessel make sure that the reagents used do not damage the FFPM seal.
- The unit may only be used for its "intended use".
- Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids etc.
- Only accessories and spare parts approved by **IKA** may be used.
- The appliance may only be opened by experts.
- Note when cleaning and servicing with an open reactor that falling parts could damage the reactor vessel. When carrying out such work, always remove the anchor agitator and the dispersal tool first!
- There may be electrostatic discharges between the medium and the dispersing/ stirring tool shaft which could pose a direct danger.
- The appliance is not suitable for manual operation.
- Set up the stand in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that the stand does not start to move.
- The agitated vessels used for stirring have to be secured. Consider on a good stability of the entire structure.
- Secure the stirring vessel against twisting.
- Glass vessels must always be secured with a clamp to prevent them spinning. When working with flasks elastic couplers must be used to prevent fracture of the glass.
- Note the operating instructions of the dispersing/stirring tool and accessories.
- Do not use the appliance without a dispersing/stirring element.
- In the event of unbalance or unusual noises, switch off the appliance immediately. Replace the dispersing element. If there is no difference after the change of the dispersing tool, return it to the dealer or the manufacturer along with a description of the fault.
- Please note that the device and the bearings can become extremely hot during use.
- Do not cover the ventilation slots on the drive in order to ensure adequate cooling of the drive.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Abrasion of the device or the rotating accessories can get into the medium you are working on.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.

Intended use

The **IKA** LR2000 system is a modular mini plant reactor system. It was designed and constructed for emulating and optimising chemical reaction processes and for mixing, dispersal and homogenising processes on a model scale.

Its working volume is 500 ml to 2000 ml.

Because the FFPM seals are highly durable the medium can be kept at a constant temperature of $-10\text{ }^{\circ}\text{C}$ to $230\text{ }^{\circ}\text{C}$. Vacuum operation is possible up to 25 mbar (see "Accessories").

By adapting the torque measuring instrument VK 600 control (see accessories) rheological material properties, flow and deformation behaviour can be evaluated.

A wide range of different **IKA** laboratory equipment gives the user a choice of applications as well as matching them optimally to specific requirements.

By linking it to the "Labworldsoft" laboratory software the user has a powerful system for measuring, controlling, regulating and completely documenting at their disposal.

Examples of Application

- Manufacture of creams, lotions, emulsions, preparation of liposome in the pharmaceutical and cosmetic fields (lipstick, sun protection creams, eyeliner etc)
- Working in solids like calcium carbonate, talcum, titanium oxide and others into liquid polymers
- Working additives and solid polymer compounds into mineral oils
- Dispersing and defibrating solids and fibres in liquids and polymers
- Manufacturing milk products and dietary foods
- Extracting enzymes from biomass

Range of use

- Laboratories

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.

Unpacking

► Unpack

- Please unpack the device carefully
- In the case of any damage a fact report must be set immediately (post, rail or forwarder)

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Setting up

The **IKA**-Laboratory reactor system has been supplied in a basic version ready for operation and can be directly set up on a level work bench.

Inserting the reactor vessel:

- Loosen clamping lever on telescopic stand
- Push stirrer drive upwards by hand (supported by a protective gas spring)
- Insert reactor vessel into instrument base
- Open container collets via crank handle
- Press stirring drive down
- Fix telescopic stand with clamping lever in this position
- Close container collets
- Close all three clamps and tighten them

The stirrer drive is firmly fitted to the stirring shaft and is a permanent part of the **IKA**-laboratory reactor. Only the reactor vessels are changed for the respective operations.


Fig. 1

Commissioning

Attention has to be paid to the operating instructions of the individual components: stirrer RE 162 P4 analog and RE 162 P7 analog respectively, and if required, the ULTRA-TURRAX T 25.

At first the filled reactor vessel is put into the container support in the instrument base. Then reactor vessel and hoses for tempering liquid are connected. For this the quick-release couplings (plug nipple) with the in- and outlet hoses of the thermostat are connected. The inlet is plugged to the bottom quick-release coupling and the outlet to the top quick-release coupling. The quick-release couplings are self-sealing on both sides so that no tempering liquid can leak out.

The **IKA**-laboratory reactor LR.2000 V is designed for operations under vacuum or atmospheric pressure. The vacuum pump is connected laterally at the left-hand side of the telescopic stand (see fig. 2).

The attained vacuum can be displayed by means of the vacuum gauge LR 2000.30 or by means of the **IKA**-vacuum controller VC 2. All further installation parts such as dispersing tool, temperature sensor, flow breaker and ground-glass joints for the reactor cover are sealed with PTFE-sleeves or O-rings.

Description

The **IKA**-laboratory reactor system has a modular construction. It consists essentially of the following construction groups (see fig. 1).

- instrument base with telescopic stand
- stirrer carrier with stirring shaft bearing
- stirrer drive
- container clamping device with reactor cover
- reactor vessel

Device base, telescopic stand

Vibration insulating plates are fitted to the bottom side of the instrument base: In the centre there are thermally insulated container supports for fixing the reactor vessel. The telescopic stand is fitted at the rear. The protective gas spring, built into the telescopic stand, allows for an easy lifting and lowering of the complete reactor. The reactor can be fixed at any given height by means of a clamping device. At the rear of the stand column there is a multiple cable channel which houses the connecting cable of the stirrer drive. The multiple cable channel can also take up the connecting cable of the ULTRA-TURRAX and the temperature measuring instrument.

Stirring motor carrier, stirring shaft bearing

The stirrer carrier is clamped to the telescopic stand column. At its lateral side a bracket is fitted for the following connections:

- vacuum connection with stop cock
- vacuum gauge or vacuum controller
- aerating valve
- vacuum hose leading to the reactor cover
- threaded bore holes for stand rods

For depositing the ULTRA-TURRAX-drive, a bracket with spigot and a stand rod is laterally fitted to the right-hand side. The stirring shaft bearing is firmly fitted into the front bearing eye of the stirrer carrier. The stirring shaft which is ball bearing supported and vacuum tight, is connected firmly to the stirrer drive by means of a coupling flange. The stirring shaft with covered ball bearings is precisely supported and sealed with radial shaft sealing rings with PTFE-sealing lips for vacuum operation. In order to keep the unavoidable abrasion of the stirring shaft caused by its running-in, away from the product, a protective ring is pushed from below onto the stirring shaft. The anchor stirrer with PTFE-wipers is also plugged from below into the stirring shaft. The rotary motion is transferred to the anchor stirrer by means of a square bar, and is axially fixed with a circlip.


Fig. 2

Stirrer drive

Please pay attention to the operating instructions of the individual component „Stirring Motor“!

The stirrer drive is a firm component of the **IKA**-Laboratory Reactor System. The following **IKA**-stirrers are at disposal for this system:

RE 162 P 4 analog and **RE 162 P7 analog** (optionally also **IKA-Visc-measuring stirrer**)

The RE 162 P4 analog is preferably used for the viscosity range 1 - 20.000 mPa s.

The RE 162 P7 analog is used for viscosities of more than 20.000 up to 100.000 mPa.s.

The coupling of the stirrer is fitted to the drive shaft and screw-connected to the coupling flange at the upper end of the stirring shaft bearing.

Plase note: The stirrer cannot be switched on while the container collets are open.

Container clamping device

The container clamping device is pivoted to the stirrer carrier. Via a crank handle at the right-hand side of the telescopic stand and via a threaded spindle, the container clamping device can be opened and closed. The collets are provided with thermal insulated PTFE-coatings. Snap closures and collets are screw-connected. The reactor vessel is centered to the reactor cover by means of the collets. After closing the three snap closure straps, the reactor vessel is drawn tightly against the reactor cover with three knurled-head screws.

Reactor cover

The reactor cover is provided with four standard openings. These are machined as female parts of ground joints (3 x HNS 29/32 and 1 x HNS 14/23). There is an additional connection for generating a vacuum. All openings are sealed with O-rings. The sealing between reactor cover and reactor vessel is also by an O-ring which is firmly fitted into the cover. ULTRA-TURRAX, thermal sensor, flow breaker and other accessories can be fitted into the free openings. For this purpose the corresponding screw-connections are removed by a guide pin and unscrewed from the cover. The guide pin is enclosed in the tool set.

The reactor cover is screw-connected to the stirring shaft bearing via a thread. The cover can be easily dismantled for cleaning or sterilizing purposes (see section cleaning and sterilization).

Stirring shaft sealing and reactor cover are sealed by an O-ring.

Use of dispersing tools S 25 KV

When using the ULTRA-TURRAX T 25, the following dispersing tools can be used:

S 25 KV - 18 G, SS 25 KV - 25 G or S 25 KV - 25 F

For fitting the dispersing tool into the reactor cover, shaft support LR 2000.40 is required.

Fitting Sequence: Insert retaining ring (2) into groove of shaft (1)

Push shaft from below into shaft support (4)

Insert retaining ring (3) into upper groove

Tighten clamping screws (5) (socket wrench in tool set)

The required retaining rings are part of the supply schedule of the shaft support.


Fig. 3

Reactor vessel

The double-walled glass vessel with a useful volume of 2000 ml is made of borosilicate glass 3.3 DIN-ISO 3585 and has a nominal width of DN 150. The reactor vessel can be heated up to a maximum temperature of 230 °C via its double-jacket. The double-jacket-connections are equipped with self-sealing quick-release-couplings. (Two locking nipples are in the accessories kit). A cork ring for putting down the reactor vessel LR 2000.1 or LR 2000.3 is included in the supply schedule.

Materials in contact with product

| | |
|-----------------|-------------------------------------|
| Reactor cover: | stainless steel 1.4571 |
| Reactor vessel: | borosilicate glass 3.3 DIN-ISO 3585 |
| Anchor stirrer: | stainless steel 1.4571 |
| Wipers: | PTFE |
| O-rings: | FKM or PTFE |
| Shaft sealing: | PTFE / 1.4571 |

Cleaning and sterilization

Cleaning:

The reactor vessel can be completely removed from the reactor so that the internal space is easily accessible for cleaning.

The anchor stirrer, axially held by a circlip, can be withdrawn from the stirring shaft, if required by a slight impact. The PTFE-wipers can also be easily drawn off the anchor stirrer thus ensuring a complete cleaning of the stirrer.

If no additional instruments such as ULTRA-TURRAX, thermal sensor or flow breaker are installed, the reactor cover can be screwed off the stirring shaft bearing for cleaning purposes.

After removing the vacuum connection, use the supplied hooked spanner for loosening the reactor cover. After removing the ground section parts, the screwed-joints for additional instruments, and after drawing off the protective ring from the stirring shaft, the reactor cover can be unscrewed from the stirring shaft bearing after the collets have been opened.

Stirring shaft and protective ring can now also be cleaned easily.

If required, all O-rings can be carefully removed with a blunt tool.

Sterilization:

- a) Chemical procedure:
By germicides (formaline, phenol, alcohol etc.) many disinfection problems can be solved. A subsequent removal of the remains of the disinfectants is important.
- b) Killing germs by moist heat:
Here perched steam with an overpressure of 2 bar at 120 °C is applied.
- c) Killing germs by hot air:
Generally a hot air sterilization is carried out at 160 °C up to 190 °C.
- d) Flammings:
This method can only be applied for individual parts without bore holes as it processes outer surfaces only.

After removing the O-ring sealings, reactor cover, anchor stirrer and other parts can be sterilized by method a), b) or c).

Maintenance

The **IKA**-laboratory system requires little maintenance. Ball bearings are sealed and permanently lubricated. The following bur-nished parts should be oiled once a month (see fig. 1 and 2):

- End plate of telescopic stand
- Clamping ring of telescopic stand
- Washers on rotary point of collets
- Threaded spindle and spindle nuts of container clamping device
- Retaining clips for ULTRA-TURRAX-drive

Spare parts


When ordering spare parts, please quote type of instrument and serial number as well as ident-number of the spare part. If you require spare parts, please ask for the corresponding service papers for determining spare parts.

Repair

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

In case of repair the device has to be cleaned and free from substances which may constitute a health hazard as well as from caustic constituents.

Wear parts

| Designation | Id.No. (FKM) | Id.No. (FFPM) |
|--------------------|-----------------|------------------|
| b O-ring 17 x 2,5 | 25 104 00 | 24 983 00 |
| c O-ring 159 x 6,0 | 25 118 00 | 24 986 00 |
| d O-ring 48 x 2,0 | 25 119 00 | 24 988 00 |
| e O-ring 37 x 2,0 | 25 112 00 | 24 984 00 |
| f O-ring 26 x 2,0 | 25 115 00 | 24 985 00 |
| g Centering ring | 25 894 00 | 25 184 00 |

| Designation | Id.No. |
|-----------------------------------|-----------|
| h Plug NS 29/32 | 19 707 00 |
| i Plug NS 14/32 | 19 710 00 |
| j Wiper | 25 053 00 |
| k Wiper | 25 052 00 |
| l Container guideway | 25 069 00 |
| a PS-SEAL high-efficiency sealing | 19 206 00 |

Fig. 4

Accessories

Accessories required for vacuum operation:

| | | Ident-No. |
|---------------------------|-------------------|-----------|
| - Vacuum gauge | LR-2000.30 - FFPM | 25 094 00 |
| - Vacuum gauge | LR-2000.31 - FKM | 25 917 00 |
| alternative: | | |
| - Vacuum controller | VC 2 | 19 806 00 |
| - Chemical diaphragm pump | MZ-2C | 19 805 00 |

Further optional accessories:

| | | |
|--|-------------------|-----------|
| - Reactor vessel glass | LR-2000.1 | 25 083 00 |
| - Reactor vessel glass with bottom drain valve | LR-2000.2 | 25 096 00 |
| - Stand support | LR-2000.53 | 25 098 00 |
| - Reactor vessel stainless steel | LR-2000.3 | 25 097 00 |
| - Stirring tool with wiper | LR-2000.10 | 25 081 00 |
| - Stirring tool | LR-2000.11 | 25 095 00 |
| | | |
| - ULTRA-TURRAX | T 25 | 16 190 00 |
| - ULTRA-TURRAX | T 25 S 5 | 16 190 05 |
| - Vacuum shaft | S 25KV-18G | 23 480 00 |
| - Vacuum shaft | S 25KV-25G | 24 669 00 |
| - Vacuum shaft | S 25KV-25F | 24 040 00 |
| - Shaft support | LR-2000.40 - FFPM | 25 092 00 |
| - Shaft support | LR-2000.41 - FKM | 25 915 00 |
| - IKATRON revolution counter | DZM 5 | 16 058 00 |
| - Measuring sensor | RLS 8 x 33 | 10 711 00 |
| - IKATRON digital temper.measuring instr. | DTM 11 | 24 350 00 |
| - Measuring sensor for installation | PT 100.5 | 25 068 00 |
| - Measuring sensor support | LR-2000.60 - FFPM | 25 093 00 |
| - Measuring sensor support | LR-2000.61 - FKM | 26 294 00 |
| | | |
| - IKA -circulating thermostat | LT 5 control | 26 066 00 |
| | | |
| - Hose - metal | LT 5.20 | 26 067 00 |
| | | |
| - Solenoid valve | VC 1.3 | 21 635 00 |
| - Pump control | VC 2.4 | 24 391 00 |
| | | |
| - Chemical diaphragm pump | MZ-2C | 19 805 00 |

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

LR 2000 P

PLEASE NOTE! For vacuum operations with the laboratory reactor LR 2000 P, the following points have to be noted and/or considered:

- The connection block at the left-hand side of the stand has to be exchanged against connection block for vacuum operation.
- The fittings in the reactor cover have to be exchanged against the fittings for vacuum operation (flow breaker, vacuum connection).

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