IKA C 48

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MANUAL INSTRUCTIONS *EN* 7

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Garantie

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Guarantee

You have purchased an original IKA machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 12 months.

For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim.

You will be liable for freight costs.

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	Page
Garantie	2
Notes on safety	7
Proper usage and description of the device	7
Unpacking	8
Scope of delivery	8
Installation	8
Placing the unit in service	9
Maintenance and repair	9
Accessories	10
Technical data	10
List of spare parts	10
Illustration of spare parts	15

Notes on safety

Keep lines and fittings through which oxygen flows free of grease.

Observe the accident prevention requirements in force for the task and the work station.

The main valve of the oxygen supply should be closed when the task and the work station.

Maintenance tasks should be performed without any excess pressure on the system.

The filling time should not be set to less than 30 seconds. This will prevent combustible substances from being whirled around in the crucible.

The valve needle (29) is held in place with safety varnish to keep it from being turned unintentionally. It should only be removed if there is no excess pressure on the system.

Proper usage and description of the device

The IKA C 48 oxygen filling station is suitable for filling decomposition vessels C1/8; C 7010, C 7011, C 7012, C 7010FZ, C 7012FZ and AOD 1.1 with oxygen.

- The 2-m pressure hose, included with delivery, can easily be laid to connect to the oxygen bottle. Therefore, the location of the C 48 filling station does not depend on the oxygen bottle being in the immediate vicinity.
- The filling pressure that is set and the pressure in the decomposition vessel can be read directly from the work station, thus preventing insufficient filling quantities because of not observing the minimum filling period.
- The filling time can be adjusted in increments.
- There is no need to be constantly occupied with the oxygen bottle because the automatic non-return valve only opens if a decomposition vessel is being filled.
- A carrier plate (screws on left-handed or right-handed) is used to prepare samples.

Unpacking

Please unpack the device carefully and make note of any damage. It is important to discover any damage resulting from transportation during unpacking instead of later. If appropriate, you should immediately record damage discovered (and whether material was sent by mail, rail or courier service).

Scope of delivery

The following items are included with delivery of the IKA C 48 oxygen filling station:

One **C 48** oxygen filling station, one holding plate, one pressure hose (2m) four O rings and one instruction manual.

Installation

Connecting the pressure hose

The pressure hose (**43**), included with delivery, is designed for a maximum pressure of 40 bar. It is 2 m long and easily be laid to connect to the oxygen bottle (smallest permissible bending radius 80 mm).


It should be connected at point „A“.

Longer plastic lines are available on request.

There is an R1/4inch threading on the C29 reducing valve (accessory) of the oxygen bottle to which the other end of the pressure hose can be connected.

An appropriate adapter is included for US pressure reducers.

Decomposition vessels

IKA decomposition vessels C1/8; C 7010, C 7011, C 7012 and AOD 1.1 are provided with an identification on the base or on the cap screw that contains information including the year of manufacturing number.

Assembly of the support plate

The support plate, included with delivery, can be screwed on right or left. The cover of the decomposition vessel can be set down on the support plate. This makes it easier to load decomposition vessels.

Placing the unit in service


The required filling pressure - usually 30 bar - can be set on the pressure reducer valve of the oxygen bottle. The adjusted pressure can also be read on the pressure gauge (11).

The automatic non-return valve (33) integrated into the oxygen filling station prevents oxygen from flowing out when this is not intended.


When the task is complete, the main valve on the oxygen supply should be closed.

The decomposition vessel should be placed on the holding plate (38 or 50). By loosening the screw (3), the holding plate (38 or 50) can be adjusted so that the filling head (32) fits exactly onto the oxygen line valve of the decomposition vessel.

All you need to do then to fill up the decomposition vessel is to swing the hand lever downwards. The increase in pressure can be observed on the pressure gauge.

The factory-set filling time can be changed by turning on the valve needle (29). Turning the safety valve counter-clockwise results in a longer filling time.

The filling time should not be less than 30 seconds. This will prevent combustible substances from being whiled around in the crucible.

The valve needle (29) is held in place with safety varnish to keep it from being turned unintentionally. It should only be turned if there is no excess pressure on the system.

Maintenance and repair

The C 48 oxygen filling station requires no maintenance to operate. It is subject only to the natural aging of the components and their statistical failure rate.


When ordering spare parts, please indicate the manufacturing number on the rating plate, the device type and the item number and description of the spare parts.

Please send in only units that have been cleaned and that are free of hazardous substances for repairs. Use only cleaning agents approved by IKA to clean IKA devices. To clean the residues below, use the corresponding cleaning agents:

Dye, coloring	Isopropanol
Construction materials	Water containing detergent
Cosmetics	Water containing detergent / isopropanol
Foodstuffs	Water containing detergent
Combustible materials	Water containing detergent

For materials that are not listed above, please inquire with IKA. Wear protective gloves to clean the units. Electrical devices must not be placed in the cleaning agent to clean them. Before using another cleaning or decontamination method than the one recommended by the manufacturer, the user must make certain by consulting with the manufacturer that the intended method will not damage or destroy the unit.

It should be expected that frequent use will result in wear on the O rings in the filling head. The filling head is screwed in place and can be unscrewed after the pressure hose is removed. The lower O ring (21), which is subject to particularly high wear and tear, can be reached without disassembling the pressure hose after the cap screw (34) has been removed.

Item 23	Seal for the non-return valve
Item 21	Seal for the decomposition vessel

Accessories

- **Holding plate** For inserting decomposition vessels C7010, C7011, C7012, C7010Fz C7012FZ and AOD 1.1
- **C 29** Reducer valve

Technical data

Permissible operating pressure:	bar	Max. 40
Pressure gauge:	bar	0 - 50
Typical filling pressure:	bar	30 steplessly adjustable
Pressure hose:		2m long
Filling speed:	sec	≥ 30 (pre-set to 30 bar)
Current filling pressure:		Can be read
Oxygen valve:		Automatic non-return valve
Dimensions W x D x H:	mm	200 x 260 x 450
Weight:	kg	6

List of spare parts

Item number	Description
3	Sink screw
7	Sealing ring
11	Pressure gauge
12	Sealing ring
13	Fitting
14	Protective sleeve
16	Plastic tube
18	Rubber foot
19	Foot (bumper)
21	O ring (2 x)
22	O ring
23	O ring, V80G
28	Distributor
29	Valve needle
32	Filling head
33	Valve pin
34	Cap screw
38	Holding plate
39	Washer nut
43	Pressure hose
50	Holding plate

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