



designed for scientists



R 1388 Propeller stirrer, 3-bladed

/// Data Sheet

Flow-efficient design. For drawing the material to be mixed from the bottom to the top. Minimum shearing forces. Used at medium to high speeds.

Stirrer Ø: 140 mm

Shaft Ø: 10 mm

Shaft length: 800 mm



designed for scientists

Technical Data

Stirrer diameter [mm]	140
Shaft diameter [mm]	10
Shaft length [mm]	800
Speed max. [rpm]	400
Material in contact with medium	AISI 316L / 1.4571
Weight [kg]	0.55

