




designed for scientists



R 1300 Dissolver stirrer

/// Data Sheet

Radial flow, for drawing the material to be mixed from the top and the bottom. High turbulence, high shearing forces. Particle reduction. Used at medium to high speeds.

Stirrer Ø: 80 mm

Shaft Ø: 8 mm

Shaft length: 350 mm



designed for scientists

Technical Data

Stirrer diameter [mm]	80
Shaft diameter [mm]	8
Shaft length [mm]	350
Speed max. [rpm]	2000
Weight [kg]	0.181

