



designed for scientists



## R 1342 Propeller stirrer, 4-bladed

/// Data Sheet

Standard stirring element. For drawing the material to be mixed from the top to the bottom. Local shearing forces. Generates axial flow in the vessel. Used at medium to high speeds.

Stirrer Ø: 50 mm

Shaft Ø: 8 mm

Shaft length: 350 mm



designed for scientists

### Technical Data

Stirrer diameter [mm]	50
Shaft diameter [mm]	8
Shaft length [mm]	350
Speed max. [rpm]	2000
Material in contact with medium	AISI 316L / 1.4571
Weight [kg]	0.142

