




designed for scientists



R 1382 Propeller stirrer, 3-bladed

/// Data Sheet

Flow-efficient design. For drawing the material to be mixed from the top and the bottom. Minimum shearing forces. Used at medium to high speeds.

Stirrer Ø: 55 mm

Shaft Ø: 8 mm

Shaft length: 350 mm



designed for scientists

Technical Data

Stirrer diameter [mm]	55
Shaft diameter [mm]	8
Shaft length [mm]	350
Speed max. [rpm]	2000
Material in contact with medium	AISI 316L / 1.4571
Weight [kg]	0.174

