

IKA

designed for scientists



H 135.601 Flask inlay 2000 ml

/// Data Sheet

Round bottom

Volume: 2000 ml

Inner diameter: 162 mm

www.ika.com

Subject to technical changes



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Technical Data

Outer diameter [mm]	195.0
Bottom shape	spherical
Useable volume max. [ml]	2000
Number of holes	1
Hole diameter [mm]	162.0
Hole depth [mm]	59.3
Hole shape [mm]	spherical
Working temperature max. [°C]	300
Material in contact with medium	Aluminum anodized

