

IKA

designed for scientists



H 135.701 Flask inlay 3000 ml

/// Data Sheet

Round bottom

Volume: 3000 ml

Inner diameter: 194.5 mm

www.ika.com

Subject to technical changes



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Technical Data

Outer diameter [mm]	235.0
Bottom shape	spherical
Useable volume max. [ml]	3000
Number of holes	1
Hole diameter [mm]	194.5
Hole depth [mm]	78.5
Hole shape [mm]	spherical
Working temperature max. [°C]	300
Material in contact with medium	Aluminum anodized