

IKA

designed for scientists



H 135.60 Flask Carrier 3000 ml without Handle

/// Data Sheet

Round bottom

Volume: 3000 ml

Max. outer diameter: 218.3 mm

www.ika.com

Subject to technical changes





designed for scientists

Technical Data

Outer diameter [mm]	218.3
Bottom shape	hollow cylindrical
Useable volume max. [ml]	3000
Number of holes	1
Hole diameter [mm]	188.0
Hole depth [mm]	66.5
Hole shape [mm]	spherical
Working temperature max. [°C]	300
Material in contact with medium	Aluminum anodized

