


IKA

designed for scientists



H 135.25 Flask carrier 250 ml without handle

/// Fiche technique

Porte-ballon 250 ml sans poignée

Pour ballons ronds

Volume : 250 ml

Diamètre extérieur : 142 mm

www.ika.com

Apporter des modifications techniques



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Données techniques

Diamètre extérieur [mm]	142.0
Bottom shape	hollow cylindrical
Volume utile max. [ml]	250
Number of holes	1
Hole diameter [mm]	88.0
Hole depth [mm]	38.5
Hole shape [mm]	spherical
Température de travail max. [°C]	300
matériaux	Aluminum anodized
Poids [kg]	0.806

