



designed for scientists

HELISTAND

/// Data Sheet

Electrical stand for ROTAVISC



designed for scientists

Technical Data

Stroke max. [mm]	142
Diameter [mm]	16
Dynamic load [kg]	3.2
Overload protection	yes
Height adjustable safety switch	yes
Dimensions (W x H x D) [mm]	351 x 343 x 252
Weight [kg]	5.5
Permissible ambient temperature [°C]	0 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 40
Voltage [V]	220 - 230 / 100 - 115
Frequency [Hz]	50/60
Power input [W]	7
Power input max. [W]	7
Power input standby [W]	0