




designed for scientists



R 1333 Anchor stirrer

/// Data Sheet

Tangential flow, high shearing rate at edges, minimum deposits on the vessel wall. Used at low speeds. Polymer reactions, even distribution of high mineral contents in liquids. The ideal stirrer for medium to highly viscous fluids. Stirrer

Ø: 150 mm

Shaft Ø: 10 mm

Shaft length: 550 mm



designed for scientists

Technical Data

Stirrer diameter [mm]	150
Shaft diameter [mm]	10
Shaft length [mm]	550
Speed max. [rpm]	800
Material in contact with medium	AISI 316L / 1.4571
Weight [kg]	0.638

