




designed for scientists



## AS 130.2 Fixing clip attachment

/// Data Sheet

For processing round flasks, measuring flasks and Erlenmeyer flasks. Please order fixing clips separately.

Dimensions (W x D x H): 230 x 230 x 24 mm

Weight: 650 g

Capacity: 20 x AS 2.1 (25 ml), 12 x AS 2.2 (50 ml), 12 x AS 2.3 (100 ml), 4 x AS 2.4 (250 ml), 4 x AS 2.5 (500 ml)

Weight: 650 g



designed for scientists

### Technical Data

Width [mm]	230
Depth [mm]	230
Weight [kg]	0.655

