

# IKA

designed for scientists



## topolino

/// Data Sheet

Extremely convenient magnetic mini-stirrer for mixing quantities up to 250 ml.

- Durable, brushless motor
- Infinitely variable speed range
- High magnetic adhesion

[www.ika.com](http://www.ika.com)

Subject to technical changes



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

### Technical Data

Number of stirring positions	1
Stirring quantity max. per stirring position (H <sub>2</sub> O) [l]	0.25
Motor rating output [W]	0.8
Speed control	stepless
Speed range [rpm]	500 - 1800
Stirring bar length max. [mm]	30
Set-up plate material	synthetic (PP)
Set-up plate dimensions [mm]	Ø 80
Dimensions (W x H x D) [mm]	95 x 37 x 115
Weight [kg]	0.484
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 21
Voltage [V]	100 - 240
Frequency [Hz]	50/60
Power input [W]	10

