



topolino mobil

/// Data Sheet

Mobile magnetic mini-stirrer for mixing quantities up to 250 ml for battery and mains operation.

- Durable, brushless motor
- Continuously adjustable speed range
- High magnetic adhesion
- Portable operation with long operating time (8-12 h)
- Short charging time (2-3 h)
- Standard replaceable AA rechargeable batteries
- Optional power mode:
 - a) Mains-free with standard batteries



designed for scientists

- b) With supplied mains adapter (without batteries)
- c) Combined mains/battery operation (with batteries fitted)





designed for scientists

Technical Data

Number of stirring positions	1
Stirring quantity max. per stirring position (H ₂ O) [l]	0.25
Speed range [rpm]	500 - 1800
Stirring bar length max. [mm]	30
Set-up plate material	synthetic (PP)
Set-up plate dimensions [mm]	Ø 80
Battery operation	yes
Dimensions (W x H x D) [mm]	140 x 42 x 140
Weight [kg]	0.747
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 20
Voltage [V]	100 - 240
Frequency [Hz]	50/60
DC Voltage [V=]	9
Current consumption [mA]	1000