




designed for scientists



IKAFLON® 50 Set (5 Pcs) slide round PTFE-coated Magnetic stirring bars

/// Data Sheet

Dimensions (L x Ø): 50 x 8 mm

Cylindrical magnetic stir bars provide outstanding centering and consistent mixing.

The integrated pivot ring reduces vibrations and wear during the stirring process.

The PTFE coating is chemically inert and resistant to aggressive chemicals

For use in containers with uneven or curved bottoms. The pivote ring maintains the optimum stirring position.



designed for scientists

Technical Data

Length [mm]	50
Diameter [mm]	8
Shape	slider round
Coating	PTFE
Material in contact with medium	Alnico V
Weight [kg]	0.05

