




designed for scientists



IKAFLON® 25 Set (5 Pcs) ellipse PTFE-coated Magnetic stirring bars

/// Data Sheet

Dimensions (L x Ø): 25 x 12 mm

Elliptical magnetic stir bars are well-suited for use in round beakers. The elliptical form ensures homogenous mixing. Minimal contact when used in plastic containers.

The PTFE coating is chemically inert and resistant to aggressive chemicals.



designed for scientists

Technical Data

Length [mm]	25
Diameter [mm]	12
Shape	ellipse
Coating	PTFE
Material in contact with medium	Alnico V
Weight [kg]	0.031

