




designed for scientists



## **IKAFLON® 25 Set (5 Pcs) cross PTFE-coated Magnetic stirring bars**

/// Data Sheet

Dimensions (L x Ø): 25 x 25 mm

Stable and quiet stirring even at high speeds.

The cross shape creates deep turbulence for fast and efficient mixing.

Even at low speeds efficient for biological applications.

Turns on the bottom of the vessel to prevent sedimentation.

The PTFE coating is chemically inert and resistant to aggressive chemicals.



designed for scientists

### Technical Data

Length [mm]	25
Diameter [mm]	25
Shape	cross
Coating	PTFE
Material in contact with medium	Alnico V
Weight [kg]	0.026

