

IKA®

Centrifuge IKA® mini G



German technology
made in the USA

IKA® mini G | Compact yet powerful!



IKA® introduces a brand new mini centrifuge ideal for applications which do not require ultra high speeds e.g., micro-filtrations and cell separation. The unit features a compact design and a small footprint which requires minimal bench space. A transparent cover allows for constant observation and it can be used for both micro-centrifuge and PCR strip tubes.

Technical data	IKA® mini G
Capacity	8 x 2.0 ml
Permitted density	1.2 kg / dm ³
Centrifugal acceleration	2000 G
Max. kinetic energy	20 Nm
Motor rating input/output	12 / 8 W
Fixed speed	6000 rpm
Reversible direction of rotation	no
Speed deviation	5%
Timer	1 – 99 min
Timer display	7 segment LED
Permissible ON time	100%
Quick stop	yes
Cover opening	automatic
Rotor exchange	without tool
Dimensions (W x D x H)	155 x 175 x 105 mm
Weight	1.4 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative moisture	80%
DC Voltage	24 V
Price	\$ 367
	Ident. No. 3958000









IKA+
Special safety features
 > The centrifuge will run only when the lid is closed
 > Lid can be opened easily at the push of a button

Features
 > Fixed speed of 6000 rpm
 > Digital timer adjustable from 1 to 99 min
 > Integrated quick stop
 > Fast startup of speed
 > Compact design and quiet operation
 > Rotor exchangeable without tools
 > Easy operation with touch keypad

3 Year warranty*

* 2+1 years after registering at www.ika.com/register

-  Adjustable digital timer
-  Quiet operation
-  Special safety features

-  Integrated quick stop
-  Rotor exchangeable without tools
-  Touch keypad for easy operation

Protection class according to DIN EN 60529: IP 30



Accessories included in delivery	
Rotor 8 x 1.5 / 2 ml	① Ident. No. 3918500
Rotor 4 x strip 8 x 0.2 ml	② Ident. No. 4531300
Adapter set 0.5 ml x 8	③ Ident. No. 4470900
Adapter set 0.4 ml x 8	④ Ident. No. 4470901
Adapter set 0.2 ml x 8	⑤ Ident. No. 4470902



Prices valid until 31st of December 2016
All prices exclusive to VAT
Subject to alteration of prices
Subject to technical changes
Indications not binding for delivery

IKA® +

Ordering made easy!

For more information about our products and to place your order, please visit:

www.ika.com

201604_Centrifuge_brochure_IWW_EN_USD

FAQ

Does IKA® offer a rotor for PCR strip tubes?

Yes, a rotor for four strips 8 x 0.2 ml is included with the IKA® mini G

What is the maximum kinetic energy for the IKA® mini G?

The maximum kinetic energy is 20 Nm

What happens when the lid of the centrifuge is opened?

The centrifuge stops immediately

What tool is needed to change the rotor?

The rotor can be changed without any tool

Is a timer included in the IKA® mini G?

Yes, the timer can be adjusted from 1 – 99 min

What is the centrifugal acceleration for the IKA® mini G?

The centrifugal acceleration for the IKA® mini G is 2000 G

Does the IKA® mini G possess a digital display?

Yes, the IKA® mini G has an integrated 7 segment LED timer

Is it possible to change the speed?

No, the speed is fixed at 6000 rpm



[IKAworldwide](#) | [#lookattheblue](#)

IKA®-Works, Inc.
2635 Northchase Pkwy SE
Wilmington, NC 28405-7419
USA

Tel. +1 800-733-3037
Tel. +1 910 452-7059
Fax +1 910 452-7693

sales@ika.net
www.ika.com



designed
to work perfectly